


COOLING PLATE BCP230



| CODE | DESCRIPTION |
|------------|----------------------|
| 40-300-203 | Cooling Plate BCP230 |



In vitro diagnostic – medical device



Manufacturer: **Bio-Optica Milano S.p.A.**

Publication date 16/03/2021
Rev. 01

Cryo module used to obtain the rapid cooling of histological samples included in paraffin.

GENERAL INFORMATION

| | | |
|------------------------|------------------------------|--|
| Dimension Features | Dimensions (W x D x H) | 410 x 605 x 405 mm |
| | Cooling surface dim. (W x D) | 370 x 350 mm |
| | Weight | 24 Kg |
| Electrical connections | Power supply | 230 V |
| | Frequency | 50 ÷ 60 Hz |
| | Power | 0,6 kW |
| | Fuses | 2 fuses of 6.3 Ampere - 5x20 mm – T6.3AH250V |
| Other Connections | Water connections | Not necessary |
| | Fumes aspiration/filtration | Not necessary |

STRUCTURAL FEATURES

- Painted sheet steel chassis.
- Large stainless steel surface that offers space for storage of up to 300 paraffin blocks, placed in vertical position.
- Provided with plexiglass transparent cover and 48 mm high edge, where is concentrated most of the cooling power, to obtain a refrigerated chamber and not only a cold supporting plane
- Cooling system by compressed cycle without CFC. Refrigerant R134A 75 grams.

TECHNICAL FEATURES

Working: Power switch button ON/OFF.
Temperature: working temperature fixed at -20° C (with top cover)

INSTALLATION

Place the instrument on a stable vibration-free laboratory table with horizontal, flat table top, as far as possible vibration-free ground. Connect the plug of the power cord to the connection socket on the rear of the instrument. Plug the power cord using the cable supplied into the wall outlet (230V~ 50/60Hz).

- Do not use extension cords or adapters and do not modify the cable supplied.



IMPORTANT: Make sure that the back grid has at least a 15 cm free space in order to allow the aeration of the cooling system and the optimal system performance. In case of installing multiple units, never place the instruments with rear aeration grids against each other. If it is not possible, leave a space of at least 60 cm between a rear grid and the other.

Applicable Standards CE marked, EN 61010: 2010, IEC 61326-1: 2012, IECEE CB test certificate
The CB Certificates notarized on BLOCKCHAIN <https://blockchain.img.it/>

Publication date 16/03/2021
Rev. 01

Bio-Optica Milano S.p.A. Via San Faustino 58 - 20134 Milano
Phone +39 02.21.27.13.1 - Fax Italia +39 02.21.53.000 - Fax Export +39 02.21.54.155
www.bio-optica.it – info@bio-optica.it – claims@bio-optica.it

Vertrieb in Deutschland :
RESOLAB GmbH
Alter Rehmer Weg 7
32574 Bad Oeynhausen
Tel. 05731-8689890
Fax 05731-8689891
Email info@resolab.de
Web www.resolab.de