

# Bio - Optica

Improving Pathology

## Bio Marking Dyes

Markers for surgical margins



**Atoxic, made of polymers of natural origin**

**They do not spread in the tissue and do not stain the reagents**



**Available in seven different colours and two sizes:**



## What are Bio-Marking Dyes

**Bio-Optica Marking Dyes**  
are special inks used to mark surgical margins.

**They have to be used when the pathologist**  
**needs to identify and distinguish one or more**  
**surgical margin at the microscope.**






## 30 ml with dropper and 240 ml bottle

### Advantages:

- Atoxic, made of polymers of natural origins
- They dry in 2-3 minutes
- Ready to use: they don't need fixing steps in other solutions
- They don't spread in the tissue
- They can be applied on fresh or fixed specimens
- They don't lose their colour in solvents during fixation and tissue processing





● **Bio Marking Dyes - 30 ml bottle with dropper**




● **PACKAGING :**

7 colours - 7 pcs	05-000-030
Blue - 8 pcs	05-014-030
Black - 8 pcs	05-015-030
Green - 8 pcs	05-016-030
Yellow - 8 pcs	05-017-030
Orange - 8 pcs	05-018-030
Red - 8 pcs	05-019-030
Purple - 8 pcs	05-020-030

● **Bio Marking Dyes - 240 ml bottle**

● **PACKAGING :**

Blue - 1 pcs	05-014-240
Black - 1 pcs	05-015-240
Green - 1 pcs	05-016-240
Yellow - 1 pcs	05-017-240
Orange - 1 pcs	05-018-240
Red - 1 pcs	05-019-240
Purple - 1 pcs	05-020-240



Vertrieb in Deutschland :  
**RESOLAB GmbH**  
 Alter Rehmer Weg 7  
 32574 Bad Oeynhausen  
 Tel. 05731-8689890  
 Fax 05731-8689891  
 Email [info@resolab.de](mailto:info@resolab.de)  
 Web [www.resolab.de](http://www.resolab.de)

**Headquarter**  
 Bio-Optica Milano S.p.a.  
 via San Faustino, 58  
 20134 Milan - ITALY -

**Telephone**  
 +39 02 2127131

**Fax**  
 +39 02 2154155

**Mail**  
[info-export@bio-optica.it](mailto:info-export@bio-optica.it)

**Web**  
[www.bio-optica.it](http://www.bio-optica.it)

**Bio - Optica**