

Istoglass

Archive for microscope slides



In the photo, from the top:

- Cover
- Module for macrosections
- Module Colorteca
- Module Istobloc
- Module Istoglass
- Base

CODE	DESCRIPTION	PACKAGING
03-5000-14	Istoglass, archive for microscope slides	1 pCS



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 22/10/2018
Rev. 02

Bio-Optica Milano S.p.A. Via San Faustino 58 - 20134 Milano
Phone +39 02.21.27.13.1 - Fax Italia +39 02.21.53.000 - Fax Export +39 02.21.54.155
www.bio-optica.it – info@bio-optica.it – claims@bio-optica.it

Vertrieb in Deutschland :
RESOLAB GmbH
Alter Rehmer Weg 7
32574 Bad Oeynhausen
Tel. 05731-8689890
Fax 05731-8689891
Email info@resolab.de
Web www.resolab.de

The filing system Istoglass is an archive realized in sheet metal, thickness 9/10, japanned white. Istoglass has fourteen drawers running on guides. Each drawer is furnished with a knob and the possibility to apply a paper label.

It is possible to insert 14 metal spring (code 03-5000-MD), one for each drawer, in order to keep standing up the slides; every Istoglass may contain up to 5.000 slides; with the separator springs the capacity is about 2.000 slides.

Each Istoglass could be equipped with a base (size 490 x 490 x 60 mm) and a cover (size 490 x 490 x 25 mm). They are realized in sheet metal, thickness 9/10, japanned white.

It is possible overlap different Istoglass elements (until maximum 10 elements). In this way, it is possible decreasing the archives' size, assuring the content protection.

Optional accessories

CODE	DESCRIPTION	PACKAGING
03-5000-BA	Base for Istoglass, white colour	1 pcs
03-5000-CO	Cover for Istoglass, white colour	1 pcs
03-5000-MD	Metal spring	1 pcs

Technical Specifications

Technical data	Material	Metal sheet
	Dimensions (l x p x h)	490 x 490 x 140 mm
	Weight (empty)	13 kg
	Weight (full)	42 kg
	Drawers per module	14
	Capacity per module	About 5000 slides
Packaging	Primary container	Cardboard box
	Lable	Paper adhesiv lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 22/10/2018

Rev. 02

Bio-Optica Milano S.p.A. Via San Faustino 58 - 20134 Milano
 Phone +39 02.21.27.13.1 - Fax Italia +39 02.21.53.000 - Fax Export +39 02.21.54.155
www.bio-optica.it – info@bio-optica.it – claims@bio-optica.it

Vertrieb in Deutschland :
 RESOLAB GmbH
 Alter Rehmer Weg 7
 32574 Bad Oeynhausen
 Tel. 05731-8689890
 Fax 05731-8689891
 Email info@resolab.de
 Web www.resolab.de

Istobloc

Archive for paraffin blocks



In the photo, from the top:

- Cover
- Module for macrosections
- Module Colorteca
- Module Istobloc
- Module Istoglass
- Base

CODE	DESCRIPTION	PACKAGING
03-B900	Istobloc, archive for paraffin blocks	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 22/10/2018
Rev. 02

Bio-Optica Milano S.p.A. Via San Faustino 58 - 20134 Milano
Phone +39 02.21.27.13.1 - Fax Italia +39 02.21.53.000 - Fax Export +39 02.21.54.155
www.bio-optica.it – info@bio-optica.it – claims@bio-optica.it

Vertrieb in Deutschland :
RESOLAB GmbH
Alter Rehmer Weg 7
32574 Bad Oeynhausen
Tel. 05731-8689890
Fax 05731-8689891
Email info@resolab.de
Web www.resolab.de

Filing system for paraffin blocks storage.

The filing system Istobloc is realized in metal sheet, thickness 9/10, japanned white.

Istobloc has fourteen drawers running on guides. Each drawer is furnished with a knob with the possibility to apply a paper label. The filing system Istobloc, with 14 drawers, can contain up to 540 embedded rings or 860 paraffin blocks.

Attention: the cassettes may be clear from paraffin to be correctly inserted in each drawer (internal drawer's dimension is 30 mm).

Each Istobloc can be equipped with a base (03-5000-BA) and a cover (03-5000-CO); realized in metal sheet, thickness 9/10, japanned white.

It is possible overlap different Istobloc elements (maximum 10). In this way, it is possible to decrease the archive's size, assuring the content protection.

Optional accessories

CODE	DESCRIPTION	PACKAGING
03-5000-BA	Base for Istoglass, white colour	1 pcs
03-5000-CO	Cover for Istoglass, white colour	1 pcs

Technical Specifications

Technical data	Material	Metal sheet
	Dimensions (l x p x h)	490 x 490 x 90 mm
	Weight (empty)	11 kg
	Weight (full)	16 kg
	Drawers per module	14
	Capacity per module	540 rings or 860 cassettes
Packaging	Primary container	Cardboard box
	Lable	Paper adhesive lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 22/10/2018

Rev. 02

Bio-Optica Milano S.p.A. Via San Faustino 58 - 20134 Milano
 Phone +39 02.21.27.13.1 - Fax Italia +39 02.21.53.000 - Fax Export +39 02.21.54.155
www.bio-optica.it – info@bio-optica.it – claims@bio-optica.it

Vertrieb in Deutschland :
 RESOLAB GmbH
 Alter Rehmer Weg 7
 32574 Bad Oeynhausen
 Tel. 05731-8689890
 Fax 05731-8689891
 Email info@resolab.de
 Web www.resolab.de